



Product name	AIA®-600II
Generic name	Fully automated random access enzyme immunoassay system
Product Overview	This diagnostic system offers the same labor-saving benefits as the AIA®-21 in a compact benchtop package. The AIA®-600II uses the same test packs and has the same comprehensive menu as the AIA®-21. This completely automated system is both an excellent stand-alone solution for small laboratories, and an ideal backup system for the AIA®-21 in larger laboratories

Typical properties / Grades

Mode of operation	True random access
Capacity	60 tests/hour
Interface	Bidirectional RS232C
Detector system	Top-to-top fluorescence detection
Assay methodology	One-step and two-step sandwich and competitive EIA