

KX-21N

Results you

can trust!



Parameters	Whole Blood Mode WBC, RBC, HGB, HCT, MCV, MCH, MCHC, PLT, LYM%, MXD%, NEUT%, LYM#, MXD#, NEUT#, RDW-SD, RDW-CV, PDW*, MPV, P-LCR* Predilute Mode WBC, RBC, PLT, HGB, HCT, MCV, MCH, MCHC
Throughput	60 samples per hour
Sample Volume	50 µl (Whole Blood Mode) 20 µl (Predilute Mode)
Data Storage	300 complete sample results with histograms
Quality Control	2 QC-programs: \bar{X} control; Levey-Jennings 6 QC-files for all SYSMEX EIGHTCHECK-3WP parameters
Sample No.	15 digits
User Interface	Built-in thermal printer Graphic printer (parallel) Line printer (serial) Host computer (RS232) Hand-held barcode reader
Multi-language software	English, French, German, Spanish, Italian, Portuguese, Japanese, Chinese
Dimensions / Weights w x h x d [inches]/[lbs]	16.5 x 19 x 14/62

*In some areas, these parameters may be used for research or investigational purpose only.

Innovative technology

- Utilizes same technology as the SYSMEX high-end systems
- Fully compatible, accurate results
- Ideal as back-up for SYSMEX full-Diff systems

Compact and fully integrated

- Extremely small footprint
- Fits easily on a laboratory bench or table
- One single, operational module

Accurate and reliable

- SYSMEX robustness: for the best possible up-time
- For “peace of mind”: results you can trust
- Sensitive flagging to support diagnosis

Easy operation and maintenance

- Requires minimal training
- Simple menus and push-button technology

Safe and secure

- Non-toxic, biodegradable reagent system
- Only two basic reagents for all parameters
- Reliable results for clinician’s and patient’s peace of mind

Network capability

- Chartable reports
- Connectable directly to a local printer

Design and specifications may be subject to change due to further product development

SYSMEX CORPORATION
1-5-1, Wakinohama-Kaigandori, Chuo-ku
Kobe 651-0073, Japan
Phone 81-78-265-0521 · Fax 81-78-265-0530
www.sysmex.co.jp

SYSMEX AMERICA, INC.
One Nelson C. White Parkway, Mundelein, IL 60060, U.S.A.
Phone 1-847-996-4500 · Fax 1-847-996-4559
www.sysmex.com/usa

SYSMEX DO BRASIL INDUSTRIA E COMERCIO LTDA.
Rua Professor Algayr Munhoz Mader
3775-Bloco C-CIC Curitiba-Parana-Brasil-CEP 81350-010
Phone 55-41-246-1314 · Fax 55-41-248-3180

SYSMEX EUROPE GMBH
Bornbarch 1, 22848 Norderstedt, Germany
Phone 49-40-527 26-0 · Fax 49-40-527 26-100
www.sysmex-europe.com

SYSMEX DEUTSCHLAND GMBH
Bornbarch 1, 22848 Norderstedt, Germany
Phone 49-40-534 10 20 · Fax 49-40-523 23 02
www.sysmex.de

SYSMEX UK LTD.
Sunrise Parkway, Linford Wood (East)
Milton Keynes, Buckinghamshire, MK14 6QF, U.K.
Phone 44-1908-669 555 · Fax 44-1908-669 409
www.sysmex.co.uk

SYSMEX SINGAPORE PTE LTD
2 Woodlands Sector 1,
#01-06 Woodlands Spectrum
Singapore 738068, Singapore
Phone 65-6221-3629 · Fax 65-6221-3687
E-mail sysmex@sysmex.com.sg
www.sysmex.com.sg

SYSMEX SHANGHAI LTD.
5th Floor, New Shanghai International Tower
360 Pudong South Road
Shanghai, 200120 China
Phone 21-6886-3300 · Fax 21-6886-3355
www.sysmex.com.cn