

VITROS® 250 Chemistry System

Stand alone or with the VITROS® 950 for standardized testing and accurate results.




The VITROS 250 is a flexible system that can function as a single routine and STAT analyzer, a dedicated system for STATs or esoterics, or as a full-menu back-up to the VITROS 950. You get the same assay menu, features and ease of use as the VITROS 950, plus on-board dilutions, with a smaller footprint. So take a closer look and you'll see that the VITROS 250 is the ideal choice for a variety of high-productivity applications in hospitals, physician offices, clinics and satellite labs.

Versatile. Accurate. Easy-to-use. Cost-effective.

Efficient and productive in multiple environments.

- Full menu backup to VITROS 950
- Small footprint
- Customized patient report for physician office laboratories
- Optional positive patient identification (PSID™) enables LIS interfacing
- Automatic urine dilutions, on-demand serum dilutions
- Up to 250 results/hour, with 40+ on-board chemistries

The system of choice for consistently reliable, high-quality results.



- Proven MicroSlide™ technology sets the standard for superior assay precision
- Minimal interference from hemolysis, lipemia and bilirubin
- Fewer repeats, redraws and dilutions than with wet chemistry systems
- Backed by the company top rated in service and support for over 17 years

The best checks and balances on any chemistry system.

- Short sample detection/clot detection/bubble detection—ensure accuracy
- Unique Smart Metering™—ensures dependable testing integrity
- Disposable tips—ensure no carryover
- Continuous process verification—ensures the right result the first time

Easy to operate—just plug in and go.

- No deionized water, no plumbing, no drains
- Load bar-coded slide cartridges, of any lot number, in one easy step—with no reagent preparation
- Accepts any size primary/secondary containers or cups
- Small sample size requirements
- User-friendly help text is available at all times
- Biosafe operations with no infectious waste

Reduce your costs. Streamline your workflow.

- Unit-dose test slides provide predictable operating costs
- Supports multiple slide lots on board at one time for better inventory control
- Low labor requirements
- QC once daily/up to 6-month calibration stability for most assays
- True random-access testing

For more information, contact your Ortho-Clinical Diagnostics Representative or call 1-800-828-6316, prompt #1.

 Ortho-Clinical Diagnostics
a Johnson & Johnson company